



8 In 1 Cavitation Body Slimming Machine , Lipo Cavitation Machine

Our Product Introduction

for more products please visit us on hifubeauty-machine.com


Basic Information

- Place of Origin: china
- Brand Name: GLM
- Certification: CE
- Model Number: LD-6
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Cavitation Body Slimming Machine
- Item: Cavitation + RF +BIO+Vacuum Theory
- Vacuum Pressure: 0-100 (KPA)
- Input Voltage: Input Voltage
- Infrared Light Wavelength: 700nm
- R.F Frequency: 5MHZ
- High Light: Lipo Cavitation Body Slimming Machine, 700nm Infrared Body Slimming Machine, Multifunctional Body Slimming Machine
- Highlight: **8 In 1 Cavitation Body Slimming Machine, Cavitation Lipo Laser Machine, Lipo Cavitation Machine**



More Images



Product Description

8 in 1 multifunctional Infrared Light RF Big Lipo Cavitation Vacuum Roller Body Slimming Beauty Machine

Technical Parameters of 8 in 1 multifunctional Infrared Light RF Big Lipo Cavitation Vacuum Roller Body Slimming Beauty Machine: The BS-SL7 combines multiple technologies for body slimming: vacuum lipo massage, infrared laser lipolysis, ultrasound cavitation, etc, which is approved to effectively reduce the appearance of cellulite, increase circulation, relieve minor muscle aches and pains, and relax muscle spasms.

The vacuum suction is like performing gymnastics on your skin and underlying tissue; suction and motorized vibrators massage your skin in a special way, thus exercising it. Gentle suction mobilizes your connective tissue and the vibrators help to increase circulation, thus facilitating exchange of fluid and clearance of fat. Infrared laser reduces skin impedance by heating skin and radio frequency energy penetrates deeply into connective tissue, and the synergistic combination of infrared laser and conducted radio frequency increases the oxygen intracellular diffusion by heating the tissue. Ultrasound cavitation can reduce unwanted fat, by passing energy deeply into human tissue, and under high-speed vibration. The heat generated will cause fat cells been resolved, thus get body slimming result.







Velashape Series --

VELASHAPE II BODY SLIMMING MACHINE

Commercial use



Mid vacuum roller

Screen

40K Cavitation

Face vacuum
suction handle

Face vacuum
suction handle

Host

Name :V9 Velashape II body slimming machine
Input voltage:220V/50HZ 110V/60HZ



- Body slimming
- Facial lift
- Wrinkle remove
- Safe and painless

FUNCTIONAL FOR WHOLE BODY TREATMENT

No Recovery 01



Physical Therapy,
Safe And Painless **02**

Any Skin Type **03**

Easy Operation **04**



PRODUCT PARAMETER



Thigh ,back.

40KHZ cavitation head

Including rf and vacuum

Use: Belly, Arms, Buttokcs



Smale vacuum
Face Suction handle

Top View



Side View



Including Vacuum,RF and Infrared light

Use: face and neck.

PRODUCT DISPLAY

Front View





Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeaauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou